


AT Tech Tip: How to Send Overlays to the IntelliKeys within Overlay Maker 3




InterACT Team

The Interdisciplinary Augmentative Communication & Technology Team


1. Open Overlay Maker 3 and select an overlay to be sent to the IntelliKeys keyboard for student use.



2. Next, single mouse-click on *Overlay Maker* in the toolbar at the top of the screen and choose the command *Send Overlay* at the bottom of the menu, then single mouse-click.



3. Once selected, the overlay will be sent to the IntelliKeys as noted by the *Sending* dialog box onscreen detailing its progress being sent in purple.



You did it! You can now use your overlay to access your onscreen activity!